

NEW

MADE IN THE UK


FlowForm® Ultra 90

FLOWFORM® ULTRA 90 IS AN EXCITING NEW MODULAR CUSHION PROVIDING EXCELLENT PRESSURE RELIEF.

Flowform® Ultra 90 consists of a specifically designed Flowform® module surrounded by a high quality viscoelastic foam top and bottom layer. This is covered with a Tuftane® liner with a high quality black polyurethane outer cover.

The outer cover is constructed of a two-way stretch Polyurethane material that is vapour-permeable yet waterproof, hard wearing and washable

Suitable for preventative and therapeutic use up to Stages 1-2 Pressure Damage.




Waterproof Tuftane® liner protects the viscoelastic foam core

Can be used any way up or round - a practical approach for community seating and wheelchair users

At the heart of the Flowform® Ultra 90 is a module filled with Flowform® liquid, offering excellent client comfort and pressure distribution

high quality viscoelastic foam on the top and bottom of the cushion encases the Flowform® module



The NEW lightweight version weighs just 2.9kg (17"x17" variant), is only 3½" high and can be used for clients up to 20 stones (127kg).

X-Sensor® pressure mapping shows excellent pressure distribution and client immersion

