

# Netti<sup>®</sup> 4U CE PLUS

COMFORT



Tilt angle of  
**-5° to +20°**



Recline angle of  
**45°**

## Netti 4U CE Plus

4 Seat width sizes  
35/40/45/50cm each with  
8cm of overall adjustment

Maximum user weight 135 kg

Tilt angle of -5° to +20°

Recline angle of 45° (90-135°)

Chassis weight 27.5-29 kg

## Current Market Leader

3 Seat width sizes  
39/44/49cm each with up to  
10cm of overall adjustment

Maximum user weight 135 kg

Tilt angle of -1° to +25°

Recline angle of 30° (90-120°)

Chassis weight 34 kg



### Additional Features of the Netti 4U CE PLUS:

- Ball and socket headrest adjustment for ideal head positioning in a range of headrest shapes. 3 options at no extra cost.
- Overall width of chair only 18cm wider than seat width.
- Tension adjustable backrest option at no additional cost.
- Netti Dynamic System elements can be added including: headrest, footrest and the Libertas Back System 15° backward responsive pneumatic spring.
- Proportional elevating leg supports.
- Rear wheel options of 24" quick release self-propel wheels (for additional cost) or 16" transit wheels (as shown).
- Back folds down for compact transportation.