


Stryker Introduces Its Latest Pressure Ulcer Prevention Solutions

SofCare®

with Three Layer Technology



Three Layer Technology

This latest 3-layer air-channeling technology creates interlocking air cells to immerse and envelop the patient providing low tissue interface pressures. The 3-layer design provides total air support preventing the patient from bottoming out because the patient is always supported by a pocket of air.


Convenient. Complete. Effective.

Stryker's SofCare family of products allows for a convenient, complete and effective way to enhance pressure redistribution and patient comfort while improving hospital outcomes.



HeelCare® Boots

- 3-layer air-channeling technology
- Adjustable straps
- Pre-inflated



Chair Cushions

- 3-layer air-channeling technology
- Provides protection while seated
- Pre-inflated



SofCare Overlay

- 3-layer air-channeling technology
- Multiple material options
- Multiple inflation options

*A health care professional must always rely on his or her own clinical judgment when deciding which treatments and procedures to use with patients. For verifying availability of Stryker products in your area please contact your Stryker representative. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker, Pressure Ulcer Management (PUM), SofCare.

Stryker UK Limited, Stryker House, Hambridge Road, Newbury, Berkshire RG14 5EG UK
+44 (0)1635 262444 www.stryker.co.uk

For more details or to place an order, please call Sumed, our exclusive distributor for mainland Britain on **01457 890980** or visit **www.sumed.co.uk**

SUMED®

Sumed International (UK) Ltd.

Integrity House, Units 1-2 Graphite Way, Hadfield, Glossop, Derbyshire SK13 1QH.

Phone 01457 890980 | Fax 01457 890990 | sales@sumed.co.uk | www.sumed.co.uk