



Uniroll[®] Mobile Commode

Taking care of your personal hygiene is easy and convenient with our wheeled Mobile Commode.

**Maximum
user weight:
120kg
(18.9 St)**



Smooth Mobility

Featuring smooth-rolling wheels, making transitions simple



Personalised Comfort

Adjustable footplates and cushioned seating provide exceptional support



Effortless Cleaning

Easily clean with removable, wipeable surfaces



Locking Castors

Features four locking castors for secure positioning and stability



Swing Away Footplates

Comes with swing away footplates for added convenience



Flip Back Armrests

Equipped with flip back armrests to facilitate patient transfers



Versatile & Reliable

The Uniroll[®] Mobile Commode is adaptable to various personal care settings. Its user-friendly design caters to your needs at home. Crafted with care from high-quality materials, the Uniroll[®] Mobile Commode chair ensures durability. Its sturdy construction provides a reliable and secure option for individuals seeking assistance in their daily bathroom routines.



sales@sumed.co.uk



www.sumedinternational.com



01457 890 980

Features

Featuring four locking castors for secure positioning, a transfer style push handle for easy manoeuvrability, and removable armrests to facilitate patient transfers. Assembly is hassle-free, requiring no tools, and the flip-up footplates add to its convenience. With dimensions of 84 × 57 × 94.5 cm, it's built to accommodate users comfortably while ensuring easy transitions and safety during daily routines.

Specification	Detail
Frame Material	Steel
Surface Coating	Yes
Upper Frame	Detachable
Toilet Seat Material	Polypropylene (PP) hard seat with Polyvinyl Chloride (PVC) soft seat
Armrest	Flip-back type, with Polyurethane (PU) armpad
Footrest	Swingaway type, height adjustable, with plastic footplate
Front Castors	5 inches diameter, with brake
Rear Castors	5 inches diameter, with brake
Additional Features	Bucket Included
Handle Height	945 mm
Armrest	760 mm
Seat Height	530 mm
Seat Width	435 mm
Overall Length	840 mm
Overall Width	570 mm
Weight Capacity	120 kg
PRODUCT CODE	UMC123

Get in Touch

Take the first step towards a more accessible and stress-free care routine. Order the Uniroll® Mobile Commode today and experience the difference it makes in your daily life. Don't let anything stand in the way of your comfort and independence and trust in the Uniroll® Mobile Commode to provide the support you need, whenever you need it.

